
API reference for labbench

Release 0.20

Dan Kuester, Shane Allman, Paul Blanchard, Yao Ma (NIST)

Apr 29, 2019

Contents

1	Submodules	1
2	labbench.backends module	3
3	labbench.core module	27
4	labbench.data module	35
5	labbench.host module	47
6	labbench.notebooks module	55
7	labbench.util module	57
8	labbench.version module	67
9	Module contents	69
	Python Module Index	71

CHAPTER 1

Submodules

labbench.backends module

class labbench.backends.**CommandLineWrapper** (*resource=None, **settings*)

Bases: *labbench.core.Device*

Virtual device representing for interacting with a command line executable. It supports threaded data logging through standard input, standard output, and standard error pipes.

On connection, the *backend* attribute is None. On a call to *execute()*, *backend* becomes a subprocess instance. When EOF is reached on the executable's stdout, the backend is assumed terminated and is reset to None.

When *execute* is called, the program runs in a subprocess. The output piped to the command line standard output is queued in a background thread. Call *read_stdout()* to retrieve (and clear) this queued stdout.

Parameters

- **arguments** (*List ()*) – list of command line arguments to pass into the executable
- **binary_path** (*Unicode ()*) – path to the file to run
- **resource** (*Unicode ()*) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation
- **timeout** (*Float ()*) – Timeout (sec) after disconnect is called before killing the process

backend = **DisconnectedBackend()**

background (**extra_arguments, **flags*)

Run the executable in the background (returning immediately while the executable continues running).

Once the background process is running, * Retrieve standard output from the executable with *self.read_stdout* * Write to standard input *self.write_stdin* * Kill the process with *self.kill* * Check whether the process is running with *self.running*

Normally, the command line arguments are determined by * appending *extra_arguments* to the global arguments in *self.settings.arguments*, and * appending pairs of [key,value] from the *flags* dictionary to the

global flags defined with command flags in local state traits in *self.settings*

Use the *self.no_state_arguments* context manager to skip these global arguments like this:

```
with self.no_state_arguments:
    self.background(...)
```

Returns None

clear()

Clear queued standard output, discarding any contents

connect()

The `connect()` method exists to comply with the `Device` object protocol. Call the `execute()` method when connected to execute the binary.

disconnect()

Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

exception_on_stderr

Use this context manager to raise exceptions if a process outputs to standard error during background execution.

foreground(*extra_arguments, **flags)

Blocking execution of the binary at the file location `self.settings.binary_path`.

Normally, the command line arguments are determined by * appending `extra_arguments` to the global arguments in `self.settings.arguments`, and * appending pairs of [key,value] from the `flags` dictionary to the

global flags defined with command flags in local state traits in `self.settings`

Use the `self.no_state_arguments` context manager to skip these global arguments like this:

```
with self.no_state_arguments:
    self.foreground(...)
```

Returns the return code of the process after its completion

kill()

no_state_arguments

Use this context manager to disable automatic use of state traits in generating argument strings.

read_stdout(wait_for=0)

Return string output queued from stdout for a process running in the background. This clears the queue.

Returns an empty string if the command line program has not been executed or is empty. Running the command line multiple times overwrites the queue.

Returns stdout

respawn

Use this context manager to respawning background execution.

running()

Return whether the executable is currently running

Returns True if running, otherwise False

class settings(device, *args, **kws)

Bases: `labbench.core.settings`

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *arguments*: List
- *arguments_min*: Int
- *binary_path*: Unicode
- *concurrency_support*: Bool
- *resource*: Unicode
- *timeout*: Float

arguments

List()

list of command line arguments to pass into the executable

arguments_min

Int()

minimum number of extra command line arguments to pass to the executable

binary_path

Unicode()

path to the file to run

concurrency_support

Bool()

whether this `Device` implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of `mylocaltraitcls` where its trait named `parameter` is redefined to have `default_value=7`.

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

timeout

Float()

Timeout (sec) after disconnect is called before killing the process

```
class state (device, *args, **kws)
    Bases: labbench.core.HasStateTraits
```

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected

Bool()

whether the Device instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

write_stdin (*text*)

```
class labbench.backends.DotNetDevice (resource=None, **settings)
    Bases: labbench.core.Device
```

Parameters **resource** (*Unicode* ()) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation

backend = **DisconnectedBackend** ()

connect ()

DotNetDevice loads the DLL; importing in the class definition tries to load a lot of DLLs on import of *ssmdevices*, which would 1) break platforms that do not support the DLL, and 2) waste memory, and 3)

disconnect ()

Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

dll_name = **None**

library = None

class settings (*device*, **args*, ***kws*)

Bases: labbench.core.HasSettingsTraits

Container for settings traits in a Device. These **settings** are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: *Bool*
- *resource*: *Unicode*

concurrency_support

Bool()

whether this *Device* implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

class state (*device*, **args*, ***kws*)

Bases: labbench.core.HasStateTraits

Container for state traits in a Device. **Getting or setting state traits** triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected

Bool()

whether the `Device` instance is connected**classmethod getter** (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using `self._device`.

One example is to send a command defined by `trait.command`.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in `self._device`. One example is to send a command defined by `trait.command`.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

class `labbench.backends.EmulatedVISADevice` (*resource=None, **settings*)Bases: `labbench.core.Device`

Act as a VISA device without dispatching any visa commands

Parameters **resource** (`Unicode()`) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation

backend = `DisconnectedBackend()`**connect** ()

Backend implementations overload this to open a backend connection to the resource.

disconnect ()

Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

generators = {<class 'labbench.core.Bool'>: <function EmulatedVISADevice.<lambda>>, <**classmethod set_backend** (*backend_name*)

backend_name can be 'py' or 'ni'

class settings (*device, *args, **kws*)Bases: `labbench.core.HasSettingsTraits`

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

`Device.settings.resource = 'insert-your-address-string-here'`

trait attributes:

- *concurrency_support*: *Bool*
- *resource*: *Unicode*

concurrency_support

Bool()

whether this `Device` implementation supports threading**classmethod define** (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

class state (*device, *args, **kws*)Bases: `labbench.backends.state`

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*
- *identity*: *Unicode*
- *options*: *Unicode*
- *status_byte*: *Unicode*

connected

Bool()

whether the `Device` instance is connected**classmethod getter** (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by `trait.command`.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

identity
Unicode()

options
Unicode()

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

status_byte
Unicode()

class labbench.backends.LabviewSocketInterface (*resource=None, **settings*)

Bases: *labbench.core.Device*

Implement the basic sockets-based control interface for labview. This implementation uses a transmit and receive socket.

State sets are implemented by simple 'command value' strings and implemented with the 'command' keyword (like VISA strings). Subclasses can therefore implement support for commands in specific labview VI the same way as in VISA commands by assigning the commands implemented in the corresponding labview VI.

The resource is the ip address of the host where the labview script is running. Use the *tx_port* and *rx_port* attributes to set the TCP/IP ports where communication is to take place.

Parameters **resource** (*Unicode()*) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation

backend = **DisconnectedBackend()**

clear()
Clear any data present in the read socket buffer.

connect()
Backend implementations overload this to open a backend connection to the resource.

delay = 1

disconnect()
Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

read (*convert_func=None*)
Receive from the rx socket until *self.rx_buffer_size* samples are received or timeout happens after *self.timeout* seconds.

Optionally, apply the conversion function to the value after it is received.

resource = '127.0.0.1'

rx_buffer_size = 1024

rx_port = 61552

class settings (*device, *args, **kws*)
Bases: *labbench.core.settings*

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: *Bool*
- *resource*: *Unicode*

concurrency_support

`Bool()`

whether this `Device` implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

resource

`Unicode()`

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

class state (*device, *args, **kws*)

Bases: `labbench.core.state`

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected

`Bool()`

whether the `Device` instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

```
timeout = 2
```

```
tx_port = 61551
```

```
write(msg)
```

Send a string over the tx socket.

```
class labbench.backends.SerialDevice (resource=None, **settings)
```

Bases: *labbench.core.Device*

A general base class for communication with serial devices. Unlike (for example) VISA instruments, there is no standardized command format like SCPI. The implementation is therefore limited to connect and disconnect, which open or close a pyserial connection object: the *link* attribute. Subclasses can read or write with the link attribute like they would any other serial instance.

A *SerialDevice* resource string is the same as the platform-dependent *port* argument to new *serial.Serial* objects.

Subclassed devices that need state descriptors will need to implement *state_get* and *state_set* methods in order to define how the state descriptors set and get operations.

Parameters

- **baud_rate** (*Int* ()) –
- **dsrdtr** (*Bool* ()) –
- **parity** (*Bytes* ()) –
- **resource** (*Unicode* ()) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation
- **rtscts** (*Bool* ()) –
- **stopbits** (*Float* ()) –
- **timeout** (*Float* ()) –
- **write_termination** (*Bytes* ()) –
- **xonxoff** (*Bool* ()) –

```
backend = DisconnectedBackend()
```

```
connect()
```

Connect to the serial device with the VISA resource string defined in *self.settings.resource*

disconnect()

Disconnect the serial instrument

classmethod from_hwid (*hwid=None, *args, **connection_params*)

Instantiate a new `SerialDevice` from a 'hwid' resource instead of a comport resource. A hwid string in windows might look something like:

```
r'PCIVEN_8086&DEV_9D3D&SUBSYS_06DC1028&REV_213&11583659&1&B3'
```

static list_ports (*hwid=None*)

List USB serial devices on the computer

Returns list of port resource information

class settings (*device, *args, **kws*)

Bases: `labbench.core.settings`

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *baud_rate*: `Int`
- *concurrency_support*: `Bool`
- *dsrdtr*: `Bool`
- *parity*: `Bytes`
- *resource*: `Unicode`
- *rtscts*: `Bool`
- *stopbits*: `Float`
- *timeout*: `Float`
- *write_termination*: `Bytes`
- *xonxoff*: `Bool`

baud_rate

`Int()`

concurrency_support

`Bool()`

whether this `Device` implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

dsrdtr

Bool()

parity

Bytes()

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

rtscts

Bool()

stopbits

Float()

timeout

Float()

write_termination

Bytes()

xonxoff

Bool()

class state (*device, *args, **kws*)

Bases: `labbench.core.HasStateTraits`

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: Bool

connected

Bool()

whether the Device instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by `trait.command`.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

class labbench.backends.SerialLoggingDevice (*resource=None, **settings*)

Bases: *labbench.backends.SerialDevice*

Manage connection, acquisition, and data retrieval on a single GPS device. The goal is to make GPS devices controllable somewhat like instruments: maintaining their own threads, and blocking during setup or stop command execution.

Listener objects must implement an attach method with one argument consisting of the queue that the device manager uses to push data from the serial port.

Parameters

- **baud_rate** (*Int* ()) –
- **data_format** (*Bytes* ()) – Data format metadata
- **dsrdtr** (*Bool* ()) –
- **max_queue_size** (*Int* ()) – Number of bytes to allocate in the data retrieval buffer
- **parity** (*Bytes* ()) –
- **poll_rate** (*Float* ()) – Data retrieval rate from the device (in seconds)
- **resource** (*Unicode* ()) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation
- **rtscts** (*Bool* ()) –
- **stop_timeout** (*Float* ()) – Delay after a call to *stop* before terminating the runloop thread
- **stopbits** (*Float* ()) –
- **timeout** (*Float* ()) –
- **write_termination** (*Bytes* ()) –
- **xonxoff** (*Bool* ()) –

backend = *DisconnectedBackend*()

clear ()

configure ()

connect ()

Connect to the serial device with the VISA resource string defined in *self.settings.resource*

disconnect ()

Disconnect the serial instrument

fetch ()

classmethod from_hwid (*hwid=None, *args, **connection_params*)

Instantiate a new *SerialDevice* from a 'hwid' resource instead of a *comport* resource. A *hwid* string in windows might look something like:

```
r'PCIVEN_8086&DEV_9D3D&SUBSYS_06DC1028&REV_213&11583659&1&B3'
```

static list_ports (*hwid=None*)

List USB serial devices on the computer

Returns list of port resource information

running ()

class settings (*device, *args, **kws*)

Bases: `labbench.backends.settings`

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *baud_rate*: *Int*
- *concurrency_support*: *Bool*
- *data_format*: *Bytes*
- *dsrdtr*: *Bool*
- *max_queue_size*: *Int*
- *parity*: *Bytes*
- *poll_rate*: *Float*
- *resource*: *Unicode*
- *rtscts*: *Bool*
- *stop_timeout*: *Float*
- *stopbits*: *Float*
- *timeout*: *Float*
- *write_termination*: *Bytes*
- *xonxoff*: *Bool*

baud_rate

Int()

concurrency_support

Bool()

whether this `Device` implementation supports threading

data_format

Bytes()

Data format metadata

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

dsrdtr

Bool()

max_queue_size

Int()

Number of bytes to allocate in the data retrieval buffer

parity

Bytes()

poll_rate

Float()

Data retrieval rate from the device (in seconds)

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

rtscts

Bool()

stop_timeout

Float()

Delay after a call to *stop* before terminating the runloop thread

stopbits

Float()

timeout

Float()

write_termination

Bytes()

xonxoff

Bool()

start ()

Start a background thread that acquires log data into a queue.

class state (*device*, **args*, ***kws*)

Bases: `labbench.core.HasStateTraits`

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected

Bool()

whether the *Device* instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

stop()

class labbench.backends.TelnetDevice (*resource=None, **settings*)

Bases: *labbench.core.Device*

A general base class for communication devices via telnet. Unlike (for example) VISA instruments, there is no standardized command format like SCPI. The implementation is therefore limited to connect and disconnect, which open or close a pyserial connection object: the *backend* attribute. Subclasses can read or write with the backend attribute like they would any other telnetlib instance.

A SerialDevice resource string is the same as the platform-dependent *port* argument to new serial.Serial objects.

Subclassed devices that need state descriptors will need to implement *state_get* and *state_set* methods in order to define how the state descriptors set and get operations.

Parameters

- **port** (*Int()*) –
- **resource** (*Unicode()*) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation
- **timeout** (*Float()*) –

backend = *DisconnectedBackend()*

connect()

Make the telnet connection to the host defined by the string in *self.settings.resource*

disconnect()

Disconnect the telnet connection

class settings (*device*, **args*, ***kws*)

Bases: `labbench.core.settings`

Container for settings traits in a Device. These **settings** are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: *Bool*
- *port*: *Int*
- *resource*: *Unicode*
- *timeout*: *Float*

concurrency_support

Bool()

whether this *Device* implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

port

Int()

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

timeout

Float()

class state (*device*, **args*, ***kws*)

Bases: `labbench.core.HasStateTraits`

Container for state traits in a Device. **Getting or setting state traits** triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected

Bool()

whether the *Device* instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

class labbench.backends.VISADevice (*resource=None, **settings*)

Bases: *labbench.core.Device*

class VISADevice (*resource, read_termination='\n', write_termination='\n'*)

VISADevice instances control VISA instruments using a pyvisa backend. Compared to direct use of pyvisa, this style of use permits use of labbench device *state* goodies for compact, readable code, as well as type checking.

For example, the following fetches the identity string from the remote instrument:

```
with VISADevice('USB0::0x2A8D::0x1E01::SG56360004::INSTR') as instr:
    print instr.state.identity
```

This is equivalent to the more pyvisa-style use as follows:

```
inst = VISADevice('USB0::0x2A8D::0x1E01::SG56360004::INSTR')
inst.connect()
print inst.query('*IDN?')
```

Use of *inst.state* makes it possible to add callbacks to support automatic state logging, or to build a UI.

Parameters **resource** (*Unicode()*) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation

backend = *DisconnectedBackend()*

connect ()

Connect to the VISA instrument defined by the VISA resource set by *self.settings.resource*. The pyvisa backend object is assigned to *self.backend*.

Returns None

Instead of calling *connect* directly, consider using *with* statements to guarantee proper disconnection if there is an error. For example, the following sets up a connected instance:

```
with VISADevice('USB0::0x2A8D::0x1E01::SG56360004::INSTR') as inst:
    print inst.state.identity
    print inst.state.status_byte
    print inst.state.options
```

would instantiate a *VISADevice* and guarantee it is disconnected either at the successful completion of the *with* block, or if there is any exception.

disconnect ()

Disconnect the VISA instrument. If you use a *with* block this is handled automatically and you do not need to call this method.

Returns None

classmethod list_resources ()**overlap_and_block (timeout=None, quiet=False)**

Append ‘;*OPC’ to each command sent during the block, then send an ‘*OPC?’ query to block until the measurement is done

This is meant to be used in *with* blocks as follows:

```
with inst.overlap_and_block():
    inst.write('long running command 1')
    inst.write('long running command 2')
```

At the end of the *with* block, a ‘*OPC?’ query is sent that blocks until the instrument signals it has finished with the commands.

Optionally, VISA timeout exceptions can be suppressed with the *quiet* argument.

Parameters

- **timeout** – Timeout, in milliseconds, or None to use *self.backend.timeout*
- **quiet** – Suppress timeout exceptions if this evaluates as True

preset ()

Convenience function to send standard SCPI ‘*RST’

query (msg, timeout=None)

Query an SCPI command to the device with pyvisa, and return a string containing the device response.

Handles debug logging, and extra *OPC commands (when in an *overlap_and_block* context), then sends a command string to the device.

Parameters *msg* (*str*) – the SCPI command to send by VISA

Returns the response to the query from the device

classmethod set_backend (backend_name)

Set the pyvisa resource manager for all VISA objects.

Parameters *str* (*backend_name*) – ‘@ni’ (the default) or ‘@py’

Returns None

class settings (*device*, **args*, ***kws*)

Bases: `labbench.core.settings`

Container for settings traits in a Device. These **settings** are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: *Bool*
- *read_termination*: *Unicode*
- *resource*: *Unicode*
- *write_termination*: *Unicode*

concurrency_support

Bool()

whether this *Device* implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

read_termination

Unicode()

termination character to indicate end of message on receive from the instrument

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

write_termination

Unicode()

termination character to indicate end of message in messages sent to the instrument

class state (*device*, **args*, ***kws*)

Bases: `labbench.core.state`

Container for state traits in a Device. **Getting or setting state traits** triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*
- *identity*: *Unicode*
- *options*: *Unicode*
- *status_byte*: *Dict*

connected

Bool()

whether the *Device* instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

identity

Unicode()

identity string reported by the instrument

options

Unicode()

options reported by the instrument

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

status_byte

Dict()

VISA status byte reported by the instrument

class suppress_timeout (**exceptions*)

Bases: *contextlib.suppress*

Context manager to suppress timeout exceptions.

Example:

```
with inst.suppress_timeout():
    inst.write('long running command 1')
    inst.write('long running command 2')
```

If the command 1 raises an exception, then command 2 will (silently) not execute.

wait()

Convenience function to send standard SCPI ‘*WAI’

write(msg)

Write an SCPI command to the device with pyvisa.

Handles debug logging, and extra *OPC commands (when in an overlap_and_block context), then sends a command string to the device.

Parameters *msg* (*str*) – the SCPI command to send by VISA

Returns None

class labbench.backends.Win32ComDevice (*resource=None, **settings*)

Bases: *labbench.core.Device*

Basic support for calling win32 COM APIs.

The python wrappers for COM drivers still basically require that threading is performed using the windows COM API, and not the python threading. Figuring this out with win32com calls within python is not for the faint of heart. Threading support is instead realized with util.ThreadSandbox, which ensures that all calls to the dispatched COM object block until the previous calls are completed from within a background thread. Set concurrency_support=True to decide whether this thread support wrapper is applied to the dispatched Win32Com object.

Parameters

- **com_object** (*Unicode()*) –
- **concurrency_support** (*Bool()*) – whether this Device implementation supports threading
- **resource** (*Unicode()*) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation

backend = DisconnectedBackend()

connect()

Connect to the win32 com object

disconnect()

Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

class settings (*device, *args, **kws*)

Bases: *labbench.core.settings*

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *com_object*: *Unicode*
- *concurrency_support*: *Bool*
- *resource*: *Unicode*

com_object

Unicode()

concurrency_support

Bool()

whether this *Device* implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

class state (*device*, **args*, ***kws*)

Bases: *labbench.core.HasStateTraits*

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected

Bool()

whether the *Device* instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

CHAPTER 3

labbench.core module

This implementation is deeply intertwined with low-level internals of traitlets and obscure details of the python object model. Consider reading the documentation closely and inheriting these objects instead of reverse-engineering this code.

```
exception labbench.core.ConnectionError
    Bases: traitlets.traitlets.TraitError

    args

    with_traceback()
        Exception.with_traceback(tb) – set self.__traceback__ to tb and return self.

exception labbench.core.DeviceNotReady
    Bases: Exception

    args

    with_traceback()
        Exception.with_traceback(tb) – set self.__traceback__ to tb and return self.

exception labbench.core.DeviceFatalError
    Bases: Exception

    args

    with_traceback()
        Exception.with_traceback(tb) – set self.__traceback__ to tb and return self.

exception labbench.core.DeviceException
    Bases: Exception

    args

    with_traceback()
        Exception.with_traceback(tb) – set self.__traceback__ to tb and return self.

exception labbench.core.DeviceConnectionLost
    Bases: Exception
```

```
    args
    with_traceback()
        Exception.with_traceback(tb) – set self.__traceback__ to tb and return self.
exception labbench.core.DeviceStateError
    Bases: traitlets.traitlets.TraitError
    args
    with_traceback()
        Exception.with_traceback(tb) – set self.__traceback__ to tb and return self.
class labbench.core.Int (default_value=traitlets.Undefined, allow_none=False, read_only=None,
                        help=None, write_only=None, cache=None, command=None, get-
                        ter=None, setter=None, remap={}, **kwargs)
    Bases: labbench.core.TraitMixin, traitlets.traitlets.CInt
    doc_attrs = ('min', 'max', 'command', 'read_only', 'write_only', 'remap', 'cache')
    info_text = 'an int'
    validate (obj, value)

class labbench.core.Float (default_value=traitlets.Undefined, allow_none=False,
                          read_only=None, help=None, write_only=None, cache=None, com-
                          mand=None, getter=None, setter=None, remap={}, **kwargs)
    Bases: labbench.core.TraitMixin, labbench.core.CFloatSteppedTraitlet
    doc_attrs = ('min', 'max', 'step', 'command', 'read_only', 'write_only', 'remap', 'cache')
    info_text = 'a float'
    validate (obj, value)

class labbench.core.Unicode (default_value=traitlets.Undefined, allow_none=False,
                           read_only=None, help=None, write_only=None, cache=None,
                           command=None, getter=None, setter=None, remap={}, **kwargs)
    Bases: labbench.core.TraitMixin, traitlets.traitlets.CUnicode
    default_value = ''
    info_text = 'a unicode string'
    validate (obj, value)

class labbench.core.Complex (default_value=traitlets.Undefined, allow_none=False,
                           read_only=None, help=None, write_only=None, cache=None,
                           command=None, getter=None, setter=None, remap={}, **kwargs)
    Bases: labbench.core.TraitMixin, traitlets.traitlets.CComplex
    info_text = 'a complex number'
    validate (obj, value)

class labbench.core.Bytes (default_value=traitlets.Undefined, allow_none=False,
                          read_only=None, help=None, write_only=None, cache=None, com-
                          mand=None, getter=None, setter=None, remap={}, **kwargs)
    Bases: labbench.core.TraitMixin, traitlets.traitlets.CBytes
    info_text = 'a bytes object'
    validate (obj, value)
```



```
class labbench.core.CaselessBytesEnum (default_value=traitlets.Undefined,
                                         allow_none=False, read_only=None, help=None,
                                         write_only=None, cache=None, command=None,
                                         getter=None, setter=None, remap={}, **kwargs)
Bases: labbench.core.TraitMixIn, labbench.core.EnumBytesTraitlet

devalidate (obj, value)

doc_attrs = ('values', 'case_sensitive', 'command', 'read_only', 'write_only', 'remap'
info_text = 'a bytes object'

validate (obj, value)

class labbench.core.Bool (default_value=traitlets.Undefined, allow_none=False, read_only=None,
                           help=None, write_only=None, cache=None, command=None, get-
                           ter=None, setter=None, remap={}, **kwargs)
Bases: labbench.core.TraitMixIn, traitlets.traitlets.CBool

default_value = False

info_text = 'a boolean'

validate (obj, value)

class labbench.core.List (default_value=traitlets.Undefined, allow_none=False, read_only=None,
                           help=None, write_only=None, cache=None, command=None, get-
                           ter=None, setter=None, remap={}, **kwargs)
Bases: labbench.core.TraitMixIn, traitlets.traitlets.List

class_init (cls, name)
    Part of the initialization which may depend on the underlying HasDescriptors class.

    It is typically overloaded for specific types.

    This method is called by MetaHasDescriptors.__init__() passing the class (cls) and name un-
    der which the descriptor has been assigned.

default_value = []

default_value_repr ()

element_error (obj, element, validator)

error (obj, value)

info ()

instance_init (obj)
    Part of the initialization which may depend on the underlying HasDescriptors instance.

    It is typically overloaded for specific types.

    This method is called by HasTraits.__new__() and in the BaseDescriptor.
    instance_init() method of descriptors holding other descriptors.

class
    alias of builtins.list

length_error (obj, value)

make_dynamic_default ()

validate (obj, value)

validate_elements (obj, value)
```

```
class labbench.core.Dict (default_value=traitlets.Undefined, allow_none=False, read_only=None,
                        help=None, write_only=None, cache=None, command=None, get-
                        ter=None, setter=None, remap={}, **kwargs)
Bases: labbench.core.TraitMixIn, traitlets.traitlets.Dict
```

```
class_init (cls, name)
```

Part of the initialization which may depend on the underlying HasDescriptors class.

It is typically overloaded for specific types.

This method is called by `MetaHasDescriptors.__init__()` passing the class (*cls*) and *name* under which the descriptor has been assigned.

```
default_value_repr ()
```

```
element_error (obj, element, validator)
```

```
error (obj, value)
```

```
info ()
```

```
instance_init (obj)
```

Part of the initialization which may depend on the underlying HasDescriptors instance.

It is typically overloaded for specific types.

This method is called by `HasTraits.__new__()` and in the `BaseDescriptor.instance_init()` method of descriptors holding other descriptors.

```
klass = None
```

```
make_dynamic_default ()
```

```
validate (obj, value)
```

```
validate_elements (obj, value)
```

```
class labbench.core.TCPAddress (default_value=traitlets.Undefined, allow_none=False,
                                read_only=None, help=None, write_only=None, cache=None,
                                command=None, getter=None, setter=None, remap={},
                                **kwargs)
```

Bases: labbench.core.TraitMixIn, traitlets.traitlets.TCPAddress

```
info_text = 'an (ip, port) tuple'
```

```
validate (obj, value)
```

```
class labbench.core.CaselessStrEnum (default_value=traitlets.Undefined, allow_none=False,
                                      read_only=None, help=None, write_only=None,
                                      cache=None, command=None, getter=None, set-
                                      ter=None, remap={}, **kwargs)
```

Bases: labbench.core.TraitMixIn, traitlets.traitlets.CaselessStrEnum

```
doc_attrs = ('values', 'command', 'read_only', 'write_only', 'remap', 'cache')
```

```
info ()
```

Returns a description of the trait.

```
validate (obj, value)
```

```
class labbench.core.Device (resource=None, **settings)
```

Bases: object

Device is the base class common to all labbench drivers. Inherit it to implement a backend, or a specialized type of driver.

Drivers that subclass *Device* get

- device connection management via context management (the *with* statement)
- test state management for easy test logging and extension to UI
- a degree automatic stylistic consistency between drivers

Parameters

- **resource** – resource identifier, with type and format determined by backend (see specific subclasses for details)
- ****local_states** – set the local state for each supplied state key and value

Note: Use *Device* by subclassing it only if you are implementing a driver that needs a new type of backend.

Several types of backends have already been implemented as part of labbench:

- *VISADevice* exposes a *pyvisa* backend for VISA Instruments
- *CommandLineWrapper* exposes a threaded pipes backend for command line tools
- *Serial* exposes a *pyserial* backend for serial port communication
- *DotNetDevice* exposes a *pythonnet* for wrapping dotnet libraries

(and others). If you are implementing a driver that uses one of these backends, inherit from the corresponding class above, not *Device*.

Parameters **resource** (`Unicode()`) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation

backend = `DisconnectedBackend()`

it is to be set in *connect* and *disconnect* by the subclass that implements the backend.

connect ()

Backend implementations overload this to open a backend connection to the resource.

disconnect ()

Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

class settings (*device*, **args*, ***kws*)

Bases: `labbench.core.HasSettingsTraits`

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: *Bool*

- *resource*: *Unicode*

concurrency_support*Bool()*

whether this *Device* implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitscls.define(parameter=7)
```

will return a deep copy of *mylocaltraitscls* where its trait named *parameter* is redefined to have *default_value=7*.

resource*Unicode()*

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

class state (*device, *args, **kws*)

Bases: *labbench.core.HasStateTraits*

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected*Bool()*

whether the *Device* instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

`labbench.core.list_devices (depth=1)`

Look for Device instances, and their names, in the calling code context (`depth == 1`) or its callers (if `depth` in `(2,3,...)`). Checks `locals()` in that context first. If no Device instances are found there, search the first argument of the first function argument, in case this is a method in a class.

exception `labbench.core.CommandNotImplementedError`

Bases: `NotImplementedError`

args

with_traceback()

Exception.with_traceback(tb) – set `self.__traceback__` to `tb` and return `self`.

labbench.data module

class labbench.data.StateAggregator

Bases: object

Aggregate state information from multiple devices. This can be the basis for automatic database logging.

get ()

Aggregate and return the current device states as configured with *observe()*.

Returns dictionary of aggregated states. Keys are str formatted

according to the *key()* (defaults to '{device name}_{state name}'). Values are the type and value of the corresponding state of the device instance.

key (device_name, state_name)

Generate a name for a state based on the names of a device and one of its remotelets

observe (devices, changes=True, always=[], never=[])

observe_settings (devices, changes=True, always=[], never=['connected'])

Configure Each time a device setting is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to *observe_settings()* replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of settings to actively update on each call to get()
- **never** – name (or iterable of multiple names) of settings to exclude from aggregated result (overrides **:param:'always'**)

observe_states (*devices*, *changes=True*, *always=[]*, *never=['connected']*)

Configure Each time a device state is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to `observe_states()` replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of states to actively update on each call to `get()`
- **never** – name (or iterable of multiple names) of states to exclude from aggregated result (overrides **param:always**)

set_device_labels (***mapping*)

Manually choose device name for a device instance.

Parameters *mapping* (*dict*) – name mapping, formatted as {device_object: 'device name'}

Returns None

```
class labbench.data.StatesToRelationalTable (path, overwrite=False,  
                                             text_relational_min=1024,  
                                             force_relational=['host_log'],  
                                             dirname_fmt='{id}' {host_time}',  
                                             nonscalar_file_type='csv', meta-  
                                             data_dirname='metadata', tar=False,  
                                             **metadata)
```

Bases: `labbench.data.StateAggregator`

Abstract base class for loggers that queue dictionaries of data before writing to disk. This extends `StateAggregator` to support

1. queuing aggregate state of devices by lists of dictionaries;
2. custom metadata in each queued aggregate state entry; and
3. custom response to non-scalar data (such as relational databasing).

Parameters

- **path** (*str*) – path to use for the master database
- **overwrite** (*bool*) – whether to overwrite the master database if it exists (otherwise, append)

append (**args*, ***kwargs*)

Add a new row of data to the list of data that awaits write to disk, *self.pending*. The row is represented as a dictionary of pairs formatted as {'column_name': 'row_value'}. These pairs come from a combination of 1) keyword arguments passed as *kwargs*, 2) a single dictionary argument, and/or 3) state traits configured automatically with *self.observe_states*.

The first pass at forming the row is the single dictionary argument

```
row = {'name1': value1, 'name2': value2, 'name3': value3}  
db.append(row)
```


The second pass is to update with values as configured with *self.observe_states*.

Keyword arguments are passed in as

```
db.append(name1=value1, name2=value2, nameN=valueN)
```

Simple “scalar” database types like numbers, booleans, and strings are added directly to the table. Non-scalar or multidimensional values are stored in a separate file (as defined in *set_path_format()*), and the path to this file is stored in the table.

The row of data is appended to list of rows pending write to disk, *self.pending*. Nothing is written to disk until *write()*.

Parameters **copy=True** (*bool*) – When *True* (the default), a deep copy o

data is used to avoid problems with overwriting references to data if *data* is reused during test. This takes some extra time; set to *False* to skip this copy operation.

Returns the dictionary representation of the row added to *self.pending*.

clear()

Remove any queued data that has been added by append.

close()

Close the file or database connection. This is an abstract base method (it be overridden by inheriting classes)

Returns None

get()

Aggregate and return the current device states as configured with *observe()*.

Returns dictionary of aggregated states. Keys are str formatted

according to the *key()* (defaults to ‘{device name}_{state name}’). Values are the type and value of the corresponding state of the device instance.

index_label = 'id'

key(device_name, state_name)

Generate a name for a state based on the names of a device and one of its remotelets

observe(devices, changes=True, always=[], never=[])

observe_settings(devices, changes=True, always=[], never=['connected'])

Configure Each time a device setting is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to *observe_settings()* replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of settings to actively update on each call to get()
- **never** – name (or iterable of multiple names) of settings to exclude from aggregated result (overrides **param:‘always’**)

observe_states (*devices*, *changes=True*, *always=[]*, *never=['connected']*)

Configure Each time a device state is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to `observe_states()` replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of states to actively update on each call to `get()`
- **never** – name (or iterable of multiple names) of states to exclude from aggregated result (overrides **param:‘always’**)

open (*path=None*)

set_device_labels (***mapping*)

Manually choose device name for a device instance.

Parameters *mapping* (*dict*) – name mapping, formatted as {device_object: ‘device name’}

Returns None

set_path_format (*format*)

Set the path name convention for relational files that is used when a table entry contains non-scalar (multidimensional) information and will need to be stored in a separate file. The entry in the aggregate states table becomes the path to the file.

Parameters **format** – a string compatible with `str.format()`, with replacement

fields defined from the keys from the current entry of results and aggregated states.

The format string follows the syntax of python’s python’s built-in `str.format()`. You may use any keys from the table to form the path. For example, consider a scenario where aggregate device states includes *inst1_frequency* of *915e6*, and *append()* has been called as *append(dut="DUT15")*. If the current aggregate state entry includes *inst1_frequency=915e6*, then the format string ‘{dut}/{inst1_frequency}’ means relative data path ‘DUT15/915e6’.

Returns None

set_relational_file_format (*format*)

Set the format to use for relational data files.

Parameters **format** (*str*) – one of ‘csv’, ‘json’, ‘feather’, or ‘pickle’

set_row_preprocessor (*func*)

Define a function that is called to modify each pending data row before it is committed to disk. It should accept a single argument, a function or other callable that accepts a single argument (the row dictionary) and returns the dictionary modified for write to disk.

setup ()

Open the file or database connection. This is an abstract base method (it be overridden by inheriting classes)

Returns None

write()

Commit any pending rows to the master database, converting non-scalar data to data files, and replacing their dictionary value with the relative path to the data file.

Returns the number of rows written

```
class labbench.data.StatesToCSV (path,          overwrite=False,      text_relational_min=1024,
                                force_relational=['host_log'],      dirname_fmt='{id}'
                                {host_time}',      nonscalar_file_type='csv',      meta-
                                data_dirname='metadata', tar=False, **metadata)
```

Bases: `labbench.data.StatesToRelationalTable`

append(*args, **kwargs)

Add a new row of data to the list of data that awaits write to disk, *self.pending*. The row is represented as a dictionary of pairs formatted as {'column_name': 'row_value'}. These pairs come from a combination of 1) keyword arguments passed as *kwargs*, 2) a single dictionary argument, and/or 3) state traits configured automatically with *self.observe_states*.

The first pass at forming the row is the single dictionary argument

```
row = {'name1': value1, 'name2': value2, 'name3': value3}
db.append(row)
```

The second pass is to update with values as configured with *self.observe_states*.

Keyword arguments are passed in as

```
db.append(name1=value1, name2=value2, nameN=valueN)
```

Simple “scalar” database types like numbers, booleans, and strings are added directly to the table. Non-scalar or multidimensional values are stored in a separate file (as defined in *set_path_format()*), and the path to this file is stored in the table.

The row of data is appended to list of rows pending write to disk, *self.pending*. Nothing is written to disk until *write()*.

Parameters **copy=True** (*bool*) – When *True* (the default), a deep copy o

data is used to avoid problems with overwriting references to data if *data* is reused during test. This takes some extra time; set to *False* to skip this copy operation.

Returns the dictionary representation of the row added to *self.pending*.

clear()

Remove any queued data that has been added by *append*.

close()

Close the file or database connection. This is an abstract base method (it be overridden by inheriting classes)

Returns None

get()

Aggregate and return the current device states as configured with *observe()*.

Returns dictionary of aggregated states. Keys are str formatted

according to the *key()* (defaults to '{device name}_{state name}'). Values are the type and value of the corresponding state of the device instance.

index_label = 'id'

key (*device_name*, *state_name*)

Generate a name for a state based on the names of a device and one of its remotelets

nonscalar_file_type = 'csv'

observe (*devices*, *changes=True*, *always=[]*, *never=[]*)

observe_settings (*devices*, *changes=True*, *always=[]*, *never=['connected']*)

Configure Each time a device setting is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to `observe_settings()` replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of settings to actively update on each call to `get()`
- **never** – name (or iterable of multiple names) of settings to exclude from aggregated result (overrides **:param:'always'**)

observe_states (*devices*, *changes=True*, *always=[]*, *never=['connected']*)

Configure Each time a device state is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to `observe_states()` replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of states to actively update on each call to `get()`
- **never** – name (or iterable of multiple names) of states to exclude from aggregated result (overrides **:param:'always'**)

open ()

Instead of calling `open` directly, consider using `with` statements to guarantee proper disconnection if there is an error. For example, the following sets up a connected instance:

```
with StatesToCSV('my.csv') as db:
    ### do the data acquisition here
    pass
```

would instantiate a `StatesToCSV` instance, and also guarantee a final attempt to write unwritten data is written, and that the file is closed when exiting the `with` block, even if there is an exception.

set_device_labels (***mapping*)

Manually choose device name for a device instance.

Parameters *mapping* (*dict*) – name mapping, formatted as {device_object: 'device name'}

Returns None

set_path_format (*format*)

Set the path name convention for relational files that is used when a table entry contains non-scalar (multidimensional) information and will need to be stored in a separate file. The entry in the aggregate states table becomes the path to the file.

Parameters **format** – a string compatible with `str.format()`, with replacement fields defined from the keys from the current entry of results and aggregated states.

The format string follows the syntax of python's python's built-in `str.format()`. You may use any keys from the table to form the path. For example, consider a scenario where aggregate device states includes `inst1_frequency` of `915e6`, and `append()` has been called as `append(dut="DUT15")`. If the current aggregate state entry includes `inst1_frequency=915e6`, then the format string `'{dut}/{inst1_frequency}'` means relative data path `'DUT15/915e6'`.

Returns None

set_relational_file_format (*format*)

Set the format to use for relational data files.

Parameters **format** (*str*) – one of 'csv', 'json', 'feather', or 'pickle'

set_row_preprocessor (*func*)

Define a function that is called to modify each pending data row before it is committed to disk. It should accept a single argument, a function or other callable that accepts a single argument (the row dictionary) and returns the dictionary modified for write to disk.

setup ()

Open the file or database connection. This is an abstract base method (it be overridden by inheriting classes)

Returns None

write ()

Commit any pending rows to the master database, converting non-scalar data to data files, and replacing their dictionary value with the relative path to the data file.

Returns the number of rows written

```
class labbench.data.StatesToSQLite (path,      overwrite=False,  text_relational_min=1024,
                                   force_relational=['host_log'],  dirname_fmt='{id}'
                                   {host_time}',    nonscalar_file_type='csv',    meta-
                                   data_dirname='metadata', tar=False, **metadata)
```

Bases: `labbench.data.StatesToRelationalTable`

append (*args, **kwargs)

Add a new row of data to the list of data that awaits write to disk, `self.pending`. The row is represented as a dictionary of pairs formatted as `{'column_name': 'row_value'}`. These pairs come from a combination of 1) keyword arguments passed as `kwargs`, 2) a single dictionary argument, and/or 3) state traits configured automatically with `self.observe_states`.

The first pass at forming the row is the single dictionary argument

```
row = {'name1': value1, 'name2': value2, 'name3': value3}
db.append(row)
```

The second pass is to update with values as configured with `self.observe_states`.

Keyword arguments are passed in as

```
db.append(name1=value1, name2=value2, nameN=valueN)
```

Simple “scalar” database types like numbers, booleans, and strings are added directly to the table. Non-scalar or multidimensional values are stored in a separate file (as defined in `set_path_format()`), and the path to this file is stored in the table.

The row of data is appended to list of rows pending write to disk, `self.pending`. Nothing is written to disk until `write()`.

Parameters `copy=True` (*bool*) – When *True* (the default), a deep copy of

`data` is used to avoid problems with overwriting references to data if `data` is reused during test. This takes some extra time; set to *False* to skip this copy operation.

Returns the dictionary representation of the row added to `self.pending`.

clear()

Remove any queued data that has been added by `append`.

close()

Close the file or database connection. This is an abstract base method (it be overridden by inheriting classes)

Returns None

get()

Aggregate and return the current device states as configured with `observe()`.

Returns dictionary of aggregated states. Keys are str formatted

according to the `key()` (defaults to ‘{device name}_{state name}’). Values are the type and value of the corresponding state of the device instance.

index_label = 'id'

key(name, attr)

The key determines the SQL column name. `df.to_sql` does not seem to support column names that include spaces

master_filename = 'master.db'

observe(devices, changes=True, always=[], never=[])

observe_settings(devices, changes=True, always=[], never=['connected'])

Configure Each time a device setting is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to `observe_settings()` replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of settings to actively update on each call to `get()`
- **never** – name (or iterable of multiple names) of settings to exclude from aggregated result (overrides **param:‘always’**)

observe_states (*devices*, *changes=True*, *always=[]*, *never=['connected']*)

Configure Each time a device state is set from python, intercept the value to include in the aggregate state.

Device may be a single device instance, or an several devices in an iterable (such as a list or tuple) to apply to each one.

Subsequent calls to `observe_states()` replace the existing list of observed states for each device.

Parameters

- **devices** – Device instance or iterable of Device instances
- **changes** (*bool*) – Whether to automatically log each time a state is set for the supplied device(s)
- **always** – name (or iterable of multiple names) of states to actively update on each call to `get()`
- **never** – name (or iterable of multiple names) of states to exclude from aggregated result (overrides **param:always**)

open ()

Instead of calling `open` directly, consider using `with` statements to guarantee proper disconnection if there is an error. For example, the following sets up a connected instance:

```
with StatesToSQLite('my.db') as db:
    ### do the data acquisition here
    pass
```

would instantiate a `StatesToCSV` instance, and also guarantee a final attempt to write unwritten data is written, and that the file is closed when exiting the `with` block, even if there is an exception.

set_device_labels (***mapping*)

Manually choose device name for a device instance.

Parameters *mapping* (*dict*) – name mapping, formatted as {device_object: 'device name'}

Returns None

set_path_format (*format*)

Set the path name convention for relational files that is used when a table entry contains non-scalar (multidimensional) information and will need to be stored in a separate file. The entry in the aggregate states table becomes the path to the file.

Parameters *format* – a string compatible with `str.format()`, with replacement fields defined from the keys from the current entry of results and aggregated states.

The format string follows the syntax of python's python's built-in `str.format()`. You may use any keys from the table to form the path. For example, consider a scenario where aggregate device states includes `inst1_frequency` of `915e6`, and `append()` has been called as `append(dut="DUT15")`. If the current aggregate state entry includes `inst1_frequency=915e6`, then the format string `'{dut}/{inst1_frequency}'` means relative data path `'DUT15/915e6'`.

Returns None

set_relational_file_format (*format*)

Set the format to use for relational data files.

Parameters *format* (*str*) – one of 'csv', 'json', 'feather', or 'pickle'

set_row_preprocessor (*func*)

Define a function that is called to modify each pending data row before it is committed to disk. It should

accept a single argument, a function or other callable that accepts a single argument (the row dictionary) and returns the dictionary modified for write to disk.

setup()

Open the file or database connection. This is an abstract base method (it be overridden by inheriting classes)

Returns None

table_name = 'master'

write()

Commit any pending rows to the master database, converting non-scalar data to data files, and replacing their dictionary value with the relative path to the data file.

Returns the number of rows written

`labbench.data.read(path_or_buf, columns=None, nrows=None, format='auto', **kws)`

Read tabular data from a file in one of various formats using pandas.

Parameters

- **path** (*str*) – path to the data file.
- **columns** – a column or iterable of multiple columns to return from the data file, or None (the default) to return all columns
- **nrows** – number of rows to read at the beginning of the table, or None (the default) to read all rows
- **format** (*str*) – data file format, one of ['pickle', 'feather', 'csv', 'json', 'csv'], or 'auto' (the default) to guess from the file extension
- **kws** – additional keyword arguments to pass to the pandas read_<ext> function matching the file extension

Returns pandas.DataFrame instance containing data read from file

`labbench.data.read_relational(path, expand_col, master_cols=None, target_cols=None,
master_nrows=None, master_format='auto',
prepend_column_name=True)`

Flatten a relational database table by loading the table located each row of *master[expand_col]*. The value of each column in this row is copied to the loaded table. The columns in the resulting table generated on each row are downselected according to *master_cols* and *target_cols*. Each of the resulting tables is concatenated and returned.

The expanded dataframe may be very large, making downselecting a practical necessity in some scenarios.

TODO: Support for a list of *expand_col*?

Parameters

- **master** (*pandas.DataFrame*) – the master database, consisting of columns containing data and columns containing paths to data files
- **expand_col** (*str*) – the column in the master database containing paths to data files that should be expanded
- **master_cols** – a column (or array-like iterable of multiple columns) listing the master columns to include in the expanded dataframe, or None (the default) pass all columns from *master*
- **target_cols** – a column (or array-like iterable of multiple columns) listing the master columns to include in the expanded dataframe, or None (the default) to pass all columns loaded from each *master[expand_col]*

- **master_path** – a string containing the full path to the master database (to help find the relational files)
- **prepend_column_name** (*bool*) – whether to prepend the name of the expanded column from the master database

Returns the expanded dataframe

`labbench.data.to_feather(data, path)`

Write a dataframe to a feather file on disk. Any index will be moved to a column, index and column name metadata will be removed, and columns names will be changed to a string.

Parameters

- **data** – dataframe to write to disk
- **path** – path to file to write

Returns None

labbench.host module

```

class labbench.host.Host (resource=None, **settings)
    Bases: labbench.core.Device

        Parameters resource (Unicode()) – the data needed to make a connection to a device. its type
            and format are determined by the subclass implementation

    backend = DisconnectedBackend()

    connect ()
        The host setup method tries to commit current changes to the tree

    disconnect ()
        Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

    get_git_browse_url
        Unicode()

    get_git_commit
        Unicode()

    get_git_remote_url
        Unicode()

    get_hostname
        Unicode()

    get_log
        Unicode()

    get_time
        Unicode()

    log_format = '%(asctime)s.%(msecs).03d %(levelname)10s %(message)s'

    metadata ()

    class settings (device, *args, **kws)
        Bases: labbench.core.HasSettingsTraits

```

Container for settings traits in a Device. These settings are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: *Bool*
- *resource*: *Unicode*

concurrency_support

Bool()

whether this Device implementation supports threading

classmethod define (**kws)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

resource

Unicode()

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

class state (device, *args, **kws)

Bases: `labbench.core.state`

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*
- *git_browse_url*: *Unicode*
- *git_commit_id*: *Unicode*
- *git_remote_url*: *Unicode*
- *hostname*: *Unicode*

- *log*: *Unicode*
- *time*: *Unicode*

connected

Bool()

whether the `Device` instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by *trait.command*.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

git_browse_url

Unicode()

git_commit_id

Unicode()

git_remote_url

Unicode()

hostname

Unicode()

log

Unicode()

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by *trait.command*.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

time

Unicode()

time_format = '%Y-%m-%d %H:%M:%S'

class labbench.host.**Email** (*resource=None, **settings*)

Bases: *labbench.core.Device*

Sends a notification message on disconnection. If an exception was thrown, this is a failure subject line with traceback information in the main body. Otherwise, the message is a success message in the subject line. Stderr is also sent.

Parameters

- **failure_message** (*Unicode()*) – subject line for test failure emails, or None to suppress success emails
- **recipients** (*List()*) –
- **resource** (*TCPAddress()*) –

- **sender** (`Unicode()`) –
- **success_message** (`Unicode()`) – subject line for test success emails, or `None` to suppress success emails

backend = `DisconnectedBackend()`

connect ()

Backend implementations overload this to open a backend connection to the resource.

disconnect ()

Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

send_summary ()

class settings (*device*, **args*, ***kws*)

Bases: `labbench.core.settings`

Container for settings traits in a Device. These **settings** are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: `Bool`
- *failure_message*: `Unicode`
- *recipients*: `List`
- *resource*: `TCPAddress`
- *sender*: `Unicode`
- *success_message*: `Unicode`

concurrency_support

`Bool()`

whether this `Device` implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value*=7.

failure_message

`Unicode()`

subject line for test failure emails, or `None` to suppress success emails

recipients

`List()`

resource
TCPAddress()

sender
Unicode()

success_message
Unicode()

subject line for test success emails, or None to suppress success emails

class state (*device, *args, **kws*)
Bases: `labbench.core.HasStateTraits`

Container for state traits in a Device. Getting or setting state traits triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:

```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected
Bool()

whether the `Device` instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by `trait.command`.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by `trait.command`.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

class `labbench.host.LogStderr` (*resource=None, **settings*)
Bases: `labbench.core.Device`

This “Device” logs a copy of messages on `sys.stderr` while connected.

Parameters **resource** (`Unicode()`) – the data needed to make a connection to a device. its type and format are determined by the subclass implementation

backend = `DisconnectedBackend()`

connect ()

Backend implementations overload this to open a backend connection to the resource.

disconnect ()

Backend implementations must overload this to disconnect an existing connection to the resource encapsulated in the object.

flush ()

log = ''

class settings (*device*, **args*, ***kws*)

Bases: `labbench.core.HasSettingsTraits`

Container for settings traits in a Device. These **settings** are stored only on the host; setting or getting these values do not trigger live updates (or any communication) with the device. These define connection addressing information, communication settings, and options that only apply to implementing python support for the device.

The device uses this container to define the keyword options supported by its `__init__` function. These are applied when you instantiate the device. After you instantiate the device, you can still change the setting with:

```
Device.settings.resource = 'insert-your-address-string-here'
```

trait attributes:

- *concurrency_support*: *Bool*
- *resource*: *Unicode*

concurrency_support

`Bool()`

whether this `Device` implementation supports threading

classmethod define (***kws*)

Return a copy of this class with default values of the traits redefined according to each keyword argument. For example:

```
mylocaltraitcls.define(parameter=7)
```

will return a deep copy of *mylocaltraitcls* where its trait named *parameter* is redefined to have *default_value=7*.

resource

`Unicode()`

the data needed to make a connection to a device. its type and format are determined by the subclass implementation

class state (*device*, **args*, ***kws*)

Bases: `labbench.core.HasStateTraits`

Container for state traits in a Device. **Getting or setting state traits** triggers live updates: communication with the device to get or set the value on the Device. Therefore, getting or setting state traits needs the device to be connected.

To set a state value inside the device, use normal python assignment:


```
device.state.parameter = value
```

To get a state value from the device, you can also use it as a normal python variable:

```
variable = device.state.parameter + 1
```

trait attributes:

- *connected*: *Bool*

connected

Bool()

whether the `Device` instance is connected

classmethod getter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The getter should take one argument: the instance of the trait to get. It should perform any operation needed to retrieve the current value of the device state corresponding to the supplied trait, using *self._device*.

One example is to send a command defined by `trait.command`.

The function should return a value that is the state from the device.

A trait that has its own getter defined will ignore this one.

classmethod setter (*func*)

Use this as a decorator to define a setter function for all traits in this class. The setter should take two arguments: the instance of the trait to get, and the value to set. It should perform any operation needed to apply the given value to the trait's state in *self._device*. One example is to send a command defined by `trait.command`.

Any return value from the function is ignored.

A trait that has its own setter defined will ignore this one.

write (*what*)

labbench.notebooks module

class labbench.notebooks.panel

Bases: object

Show tables summarizing device settings and states in jupyter notebook. Only a single panel will be shown in a python kernel.

Parameters **source** – Either an integer indicating how far up the calling tree to search

for Device instances, or a *labbench.Testbed* instance. :param ncols: Maximum number of devices to show on each row

children = []

devices = {}

ncols = 2

widget = None

labbench.notebooks.log_progress(*sequence*, *every=None*, *size=None*, *title=None*)

Indicate slow progress through a long sequence.

This code is adapted here from <https://github.com/alexanderkuk/log-progress> where it was provided under the MIT license.

Parameters

- **sequence** – iterable to monitor
- **every** – the number of iterations to skip between updating the progress bar, or None to update all
- **size** – number of elements in the sequence (required only for generators with no length estimate)
- **title** – title text

Returns iterator that yields the elements of *sequence*

`labbench.notebooks.range(*args, **kws)`

the same as python *range*, but with a progress bar representing progress iterating through the range

`labbench.notebooks.linspace(*args, **kws)`

the same as `numpy.linspace`, but with a progress bar representing progress iterating through the range, and an optional `title=` keyword argument to set the title

labbench.util module

labbench.util.concurrently(*objs, **kws)

If **objs* are callable (like functions), call each of **objs* in concurrent threads. If **objs* are context managers (such as Device instances to be connected), enter each context in concurrent threads.

Multiple references to the same function in *funcs* only result in one call. The *catch* and *flatten* arguments may be callables, in which case they are executed (and each flag value is treated as defaults).

Parameters

- **objs** – each argument may be a callable (function or class that defines a `__call__` method), or context manager (such as a Device instance)
- **catch** – if *False* (the default), a *ConcurrentException* is raised if any of *funcs* raise an exception; otherwise, any remaining successful calls are returned as normal
- **flatten** – if not callable and evaluates as *True*, updates the returned dictionary with the dictionary (instead of a nested dictionary)
- **nones** – if not callable and evaluates as *True*, includes entries for calls that return *None* (default is *False*)
- **traceback_delay** – if *False*, immediately show traceback information on a thread exception; if *True* (the default), wait until all threads finish

Returns the values returned by each function

Return type dictionary of keyed by function

Here are some examples:

Example Call each function *myfunc1* and *myfunc2*, each with no arguments:

```
>>> def do_something_1 ():
>>>     time.sleep(0.5)
>>>     return 1
>>> def do_something_2 ():
>>>     time.sleep(1)
```

(continues on next page)

(continued from previous page)

```

>>>     return 2
>>> rets = concurrent(myfunc1, myfunc2)
>>> rets[do_something_1]
1

```

Example To pass arguments, use the Call wrapper

```

>>> def do_something_3 (a,b,c):
>>>     time.sleep(2)
>>>     return a,b,c
>>> rets = concurrent(myfunc1, Call(myfunc3,a,b,c=c))
>>> rets[do_something_3]
a, b, c

```

Caveats

- Because the calls are in different threads, not different processes, this should be used for IO-bound functions (not CPU-intensive functions).
- Be careful about thread safety.

When the callable object is a Device method, `:func concurrency:` checks the Device object `state.concurrency_support` for compatibility before execution. If this check returns *False*, this method raises a `ConcurrentException`.

`labbench.util.sequentially(*funcs, **kws)`

Call each function or method listed in **funcs* sequentially. The goal is to emulate the behavior of the *concurrently* function, with some of the same support for updating result dictionaries.

Multiple references to the same function in *funcs* only result in one call. The *catch* and *flatten* arguments may be callables, in which case they are executed (and their values are treated as defaults).

Parameters *objs* – each argument may be a callable (function or class that defines a `__call__` method), or context manager (such as a Device instance) :param *catch*: if *False* (the default), a `ConcurrentException` is raised if any of *funcs* raise an exception; otherwise, any remaining successful calls are returned as normal :param *flatten*: if not callable and evaluates as *True*, updates the returned dictionary with the dictionary (instead of a nested dictionary) :param *nones*: if not callable and evaluates as *True*, includes entries for calls that return *None* (default is *False*) :return: the values returned by each function :rtype: dictionary of keyed by function.

Here are some examples:

Example Call each function *myfunc1* and *myfunc2*, each with no arguments:

```

>>> import labbench as lb
>>> def do_something_1 ():
>>>     time.sleep(0.5)
>>>     return 1
>>> def do_something_2 ():
>>>     time.sleep(1)
>>>     return 2
>>> rets = lb.sequentially(myfunc1, myfunc2)
>>> rets[do_something_1]
1

```

Example To pass arguments, use the Call wrapper

```

>>> def do_something_3 (a,b,c):
>>>     time.sleep(2)
>>>     return a,b,c
>>> rets = lb.sequentially(myfunc1, Call(myfunc3,a,b,c=c))
>>> rets[do_something_3]
a, b, c

```

Because `:func sequentially:` does not use threading, it does not check whether a Device method supports concurrency before it runs.

class `labbench.util.Call` (*func, *args, **kws*)

Bases: `object`

Wrap a function to apply arguments for threaded calls to *concurrently*. This can be passed in directly by a user in order to provide arguments; otherwise, it will automatically be wrapped inside *concurrently* to keep track of some call metadata during execution.

static cleanup (*func_in*)

Cleanup threading (concurrent execution only)

set_queue (*queue*)

static setup (*func_in*)

Setup threading (concurrent execution only), including checks for whether a Device instance indicates it supports concurrent execution or not.

exception `labbench.util.ConcurrentException`

Bases: `Exception`

args

with_traceback ()

`Exception.with_traceback(tb)` – set `self.__traceback__` to `tb` and return `self`.

class `labbench.util.ConfigStore`

Bases: `object`

Define dictionaries of configuration settings in subclasses of this object. Each dictionary should be an attribute of the subclass. The `all()` class method returns a flattened dictionary consisting of all values of these dictionary attributes, keyed according to `'{attr_name}_{attr_key}'`, where `{attr_name}` is the name of the dictionary attribute and `{attr_key}` is the nested dictionary key.

classmethod `all()`

classmethod `frame()`

class `labbench.util.ConcurrentRunner`

Bases: `object`

Concurrently runs all staticmethods or classmethods defined in the subclass.

class `labbench.util.FilenameDict` (**args, **kws*)

Bases: `sortedcontainers.sorteddict.SortedDict`

Sometimes instrument configuration file can be defined according to a combination of several test parameters.

This class provides a way of mapping these parameters to and from a filename string.

They keys are sorted alphabetically, just as in the underlying `SortedDict`.

clear ()

Remove all items from sorted dict.

Runtime complexity: $O(n)$

copy()

Return a shallow copy of the sorted dict.

Runtime complexity: $O(n)$

Returns new sorted dict

classmethod from_filename(filename)

Convert from a FilenameDict filename string to a FilenameDict object.

classmethod from_index(df, value=None)

Make a FilenameDict where the keys are taken from df.index and the values are constant values provided.

classmethod fromkeys(iterable, value=None)

Return a new sorted dict initialized from *iterable* and *value*.

Items in the sorted dict have keys from *iterable* and values equal to *value*.

Runtime complexity: $O(n \cdot \log(n))$

Returns new sorted dict

get()

Return the value for key if key is in the dictionary, else default.

iloc

Cached reference of sorted keys view.

Deprecated in version 2 of Sorted Containers. Use `SortedDict.keys()` instead.

items()

Return new sorted items view of the sorted dict's items.

See `SortedItemsView` for details.

Returns new sorted items view

key

Function used to extract comparison key from keys.

Sorted dict compares keys directly when the key function is none.

keys()

Return new sorted keys view of the sorted dict's keys.

See `SortedKeysView` for details.

Returns new sorted keys view

peekitem(index=-1)

Return (key, value) pair at *index* in sorted dict.

Optional argument *index* defaults to -1, the last item in the sorted dict. Specify *index*=0 for the first item in the sorted dict.

Unlike `SortedDict.popitem()`, the sorted dict is not modified.

If the *index* is out of range, raises `IndexError`.

Runtime complexity: $O(\log(n))$

```
>>> sd = SortedDict({'a': 1, 'b': 2, 'c': 3})
>>> sd.peekitem()
('c', 3)
>>> sd.peekitem(0)
('a', 1)
```

(continues on next page)

(continued from previous page)

```
>>> sd.peekitem(100)
Traceback (most recent call last):
...
IndexError: list index out of range
```

Parameters `index` (*int*) – index of item (default -1)

Returns key and value pair

Raises **IndexError** – if *index* out of range

pop (*key*, *default*=<not-given>)

Remove and return value for item identified by *key*.

If the *key* is not found then return *default* if given. If *default* is not given then raise `KeyError`.

Runtime complexity: $O(\log(n))$ – approximate.

```
>>> sd = SortedDict({'a': 1, 'b': 2, 'c': 3})
>>> sd.pop('c')
3
>>> sd.pop('z', 26)
26
>>> sd.pop('y')
Traceback (most recent call last):
...
KeyError: 'y'
```

Parameters

- **key** – *key* for item
- **default** – *default* value if key not found (optional)

Returns value for item

Raises **KeyError** – if *key* not found and *default* not given

popitem (*index*=-1)

Remove and return (*key*, *value*) pair at *index* from sorted dict.

Optional argument *index* defaults to -1, the last item in the sorted dict. Specify *index*=0 for the first item in the sorted dict.

If the sorted dict is empty, raises `KeyError`.

If the *index* is out of range, raises `IndexError`.

Runtime complexity: $O(\log(n))$

```
>>> sd = SortedDict({'a': 1, 'b': 2, 'c': 3})
>>> sd.popitem()
('c', 3)
>>> sd.popitem(0)
('a', 1)
>>> sd.popitem(100)
Traceback (most recent call last):
...
IndexError: list index out of range
```

Parameters `index` (*int*) – index of item (default -1)

Returns key and value pair

Raises

- **KeyError** – if sorted dict is empty
- **IndexError** – if *index* out of range

setdefault (*key*, *default=None*)

Return value for item identified by *key* in sorted dict.

If *key* is in the sorted dict then return its value. If *key* is not in the sorted dict then insert *key* with value *default* and return *default*.

Optional argument *default* defaults to none.

Runtime complexity: $O(\log(n))$ – approximate.

```
>>> sd = SortedDict()
>>> sd.setdefault('a', 1)
1
>>> sd.setdefault('a', 10)
1
>>> sd
SortedDict({'a': 1})
```

Parameters

- **key** – key for item
- **default** – value for item (default None)

Returns value for item identified by *key*

update (**args*, ***kwargs*)

Update sorted dict with items from *args* and *kwargs*.

Overwrites existing items.

Optional arguments *args* and *kwargs* may be a mapping, an iterable of pairs or keyword arguments. See `SortedDict.__init__()` for details.

Parameters

- **args** – mapping or iterable of pairs
- **kwargs** – keyword arguments mapping

values ()

Return new sorted values view of the sorted dict's values.

See `SortedValuesView` for details.

Returns new sorted values view

`labbench.util.hash_caller` (*call_depth=1*)

Use introspection to return an SHA224 hex digest of the caller, which is almost certainly unique to the combination of the caller source code and the arguments passed it.

`labbench.util.kill_by_name` (**names*)

Kill one or more running processes by the name(s) of matching binaries.

Parameters **names** (*str*) – list of names of processes to kill

Example

```
>>> # Kill any binaries called 'notepad.exe' or 'notepad2.exe'
>>> kill_by_name('notepad.exe', 'notepad2.exe')
```

Notes

Looks for a case-insensitive match against the `Process.name()` in the `psutil` library. Though `psutil` is cross-platform, the naming convention returned by `name()` is platform-dependent. In windows, for example, `name()` usually ends in `‘.exe’`.

`labbench.util.check_master()`

Raise `ThreadEndedByMaster` if the master thread as requested this thread to end.

`labbench.util.retry(exception_or_exceptions, tries=4, delay=0, backoff=0, exception_func=<function <lambda>>)`

This decorator causes the function call to repeat, suppressing specified exception(s), until a maximum number of retries has been attempted. - If the function raises the exception the specified number of times, the underlying exception is raised. - Otherwise, return the result of the function call.

Example

The following retries the telnet connection 5 times on `ConnectionRefusedError`:

```
import telnetlib

# Retry a telnet connection 5 times if the telnet library raises_
↳ ConnectionRefusedError
@retry(ConnectionRefusedError, tries=5)
def connect(host, port):
    t = telnetlib.Telnet()
    t.open(host, port, 5)
    return t
```

Inspired by <https://github.com/saltcrane/retry-decorator> which is released under the BSD license.

Parameters

- **exception_or_exceptions** – Exception (sub)class (or tuple of exception classes) to watch for
- **tries** (*int*) – number of times to try before giving up
- **delay** (*float*) – initial delay between retries in seconds
- **backoff** (*float*) – backoff to multiply to the delay for each retry
- **exception_func** (*callable*) – function to call on exception before the next retry

`labbench.util.show_messages(minimum_level)`

Configure screen debug message output for any messages as least as important as indicated by *level*.

Parameters `minimum_level` – One of `‘debug’`, `‘warning’`, `‘error’`, or `None`. If `None`, there will be no output.

Returns `None`

`labbench.util.sleep(seconds, tick=1.0)`

Drop-in replacement for `time.sleep` that raises `ConcurrentException` if another thread requests that all threads stop.

`labbench.util.stopwatch(desc="")`

Time a block of code using a `with` statement like this:

```
>>> with stopwatch('sleep statement') :
>>>     time.sleep(2)
sleep statement time elapsed 1.999s.
```

Parameters `desc (str)` – text for display that describes the event being timed

Returns context manager

class `labbench.util.Testbed (config=None, concurrent=True)`

Bases: object

Base class for testbeds that contain multiple Device instances or other objects like database managers that implement context management.

Use a `with` block with the testbed instance to connect everything at once like so:

```
with Testbed() as testbed:
    # use the testbed here
    pass
```

or optionally connect only a subset of devices like this:

```
testbed = Testbed()
with testbed.dev1, testbed.dev2:
    # use the testbed.dev1 and testbed.dev2 here
    pass
```

Make your own subclass of `Testbed` with a custom `make` method to define the Device or database manager instances, and a custom `startup` method to implement custom code to set up the testbed after each Device is connected.

cleanup()

This is called automatically immediately before disconnect if the testbed is connected using the `with` statement block.

Implement any custom code here in `Testbed` subclasses to implement startup of the testbed given connected Device instances.

make()

Implement this method in a subclass of `Testbed`. It should the drivers as attributes of the `Testbed` instance, for example:

```
self.dev1 = MyDevice()
```

This is called automatically when when the testbed class is instantiated.

startup()

This is called automatically after connect if the testbed is connected using the `with` statement block.

Implement any custom code here in `Testbed` subclasses to implement startup of the testbed given connected Device instances.

class `labbench.util.ThreadSandbox (factory, should_sandbox_func=None)`

Bases: object

exception `labbench.util.ThreadEndedByMaster`

Bases: `RuntimeError`

args**with_traceback()**

Exception.with_traceback(tb) – set self.__traceback__ to tb and return self.

labbench.util.until_timeout(exception_or_exceptions, timeout, delay=0, backoff=0, exception_func=<function <lambda>>)

This decorator causes the function call to repeat, suppressing specified exception(s), until the specified timeout period has expired. - If the timeout expires, the underlying exception is raised. - Otherwise, return the result of the function call.

Inspired by <https://github.com/saltycrane/retry-decorator> which is released under the BSD license.

Example

The following retries the telnet connection for 5 seconds on ConnectionRefusedError:

```
import telnetlib

@until_timeout(ConnectionRefusedError, 5)
def connect(host, port):
    t = telnetlib.Telnet()
    t.open(host, port, 5)
    return t
```

Parameters

- **exception_or_exceptions** – Exception (sub)class (or tuple of exception classes) to watch for
- **timeout** (*float*) – time in seconds to continue calling the decorated function while suppressing exception_or_exceptions
- **delay** (*float*) – initial delay between retries in seconds
- **backoff** (*float*) – backoff to multiply to the delay for each retry
- **exception_func** (*callable*) – function to call on exception before the next retry

CHAPTER 8

labbench.version module

CHAPTER 9

Module contents

I

- labbench, [69](#)
- labbench.backends, [3](#)
- labbench.core, [27](#)
- labbench.data, [35](#)
- labbench.host, [47](#)
- labbench.notebooks, [55](#)
- labbench.util, [57](#)
- labbench.version, [67](#)

-
- all()labbench.util.ConfigStore class method, 59
 - append()labbench.data.StatesToCSV method, 39
 - append()labbench.data.StatesToRelationalTable method, 36
 - append()labbench.data.StatesToSQLite method, 41
 - argslabbench.core.CommandNotImplementedError attribute, 33
 - argslabbench.core.ConnectionError attribute, 27
 - argslabbench.core.DeviceConnectionLost attribute, 27
 - argslabbench.core.DeviceException attribute, 27
 - argslabbench.core.DeviceFatalError attribute, 27
 - argslabbench.core.DeviceNotReady attribute, 27
 - argslabbench.core.DeviceStateError attribute, 28
 - argslabbench.util.ConcurrentException attribute, 59
 - argslabbench.util.ThreadEndedByMaster attribute, 64
 - argumentslabbench.backends.CommandLineWrapper.settings attribute, 5
 - arguments_minlabbench.backends.CommandLineWrapper.settings attribute, 5
 - backendlabbench.backends.CommandLineWrapper attribute, 3
 - backendlabbench.backends.DotNetDevice attribute, 6
 - backendlabbench.backends.EmulatedVISADevice attribute, 8
 - backendlabbench.backends.LabviewSocketInterface attribute, 10
 - backendlabbench.backends.SerialDevice attribute, 12
 - backendlabbench.backends.SerialLoggingDevice attribute, 15
 - backendlabbench.backends.TelnetDevice attribute, 18
 - backendlabbench.backends.VISADevice attribute, 20
 - backendlabbench.backends.Win32ComDevice attribute, 24
 - backendlabbench.core.Device attribute, 31
 - backendlabbench.host.Email attribute, 50
 - backendlabbench.host.Host attribute, 47
 - backendlabbench.host.LogStderr attribute, 51
 - background()labbench.backends.CommandLineWrapper method, 3
 - baud_ratelabbench.backends.SerialDevice.settings attribute, 13
 - baud_ratelabbench.backends.SerialLoggingDevice.settings attribute, 16
 - binary_pathlabbench.backends.CommandLineWrapper.settings attribute, 5
 - Boolclass in labbench.core, 29
 - Bytesclass in labbench.core, 28
 - Callclass in labbench.util, 59
 - CaselessBytesEnumclass in labbench.core, 28
 - CaselessStrEnumclass in labbench.core, 30
 - check_master()in module labbench.util, 63
 - childrenlabbench.notebooks.panel attribute, 55
 - class_init()labbench.core.Dict method, 30
 - class_init()labbench.core.List method, 29
 - cleanup()labbench.util.Call static method, 59
 - cleanup()labbench.util.Testbed method, 64
 - clear()labbench.backends.CommandLineWrapper method, 4
 - clear()labbench.backends.LabviewSocketInterface method, 10
 - clear()labbench.backends.SerialLoggingDevice method, 15
 - clear()labbench.data.StatesToCSV method, 39
 - clear()labbench.data.StatesToRelationalTable method, 37
 - clear()labbench.data.StatesToSQLite method, 42
 - clear()labbench.util.FilenameDict method, 59
 - close()labbench.data.StatesToCSV method, 39
 - close()labbench.data.StatesToRelationalTable method, 37
 - close()labbench.data.StatesToSQLite method, 42
 - com_objectlabbench.backends.Win32ComDevice.settings attribute, 25
 - CommandLineWrapperclass in labbench.backends, 3
 - CommandLineWrapper.settingsclass in labbench.backends, 4
 - CommandLineWrapper.stateclass in labbench.backends, 5
 - CommandNotImplementedError, 33
 - Complexclass in labbench.core, 28

[concurrency_supportlabbench.backends.CommandLineWrapper.settingslabbench.backends.CommandLineWrapper.state attribute, 5](#)
[concurrency_supportlabbench.backends.DotNetDevice.settingsconnectedlabbench.backends.DotNetDevice.state attribute, 7](#)
[concurrency_supportlabbench.backends.EmulatedVISADevice.settingslabbench.backends.EmulatedVISADevice.state attribute, 9](#)
[concurrency_supportlabbench.backends.LabviewSocketInterface.settingslabbench.backends.LabviewSocketInterface.state attribute, 11](#)
[concurrency_supportlabbench.backends.SerialDevice.settingslabbench.backends.SerialDevice.state attribute, 14](#)
[concurrency_supportlabbench.backends.SerialLoggingDevice.settingslabbench.backends.SerialLoggingDevice.state attribute, 18](#)
[concurrency_supportlabbench.backends.SerialDevice.settingsconnectedlabbench.backends.TelnetDevice.state attribute, 13](#)
[concurrency_supportlabbench.backends.SerialLoggingDevice.settingslabbench.backends.VISADevice.state attribute, 16](#)
[concurrency_supportlabbench.backends.TelnetDevice.settingsconnectedlabbench.backends.Win32ComDevice.state attribute, 19](#)
[concurrency_supportlabbench.backends.VISADevice.settingsconnectedlabbench.core.Device.state attribute, 22](#)
[concurrency_supportlabbench.backends.Win32ComDevice.settingsconnectedlabbench.host.Email.state attribute, 25](#)
[concurrency_supportlabbench.backends.Win32ComDevice.settingsconnectedlabbench.host.Host.state attribute, 25](#)
[concurrency_supportlabbench.backends.Win32ComDevice.settingsconnectedlabbench.host.LogStderr.state attribute, 53](#)
[concurrency_supportlabbench.core.Device.settings ConnectionError, 27](#)
[concurrency_supportlabbench.core.Device.settings copy\(\)labbench.util.FilenameDict method, 59](#)
[concurrency_supportlabbench.host.Email.settings attribute, 50](#)
[concurrency_supportlabbench.host.Host.settings attribute, 48](#)
[concurrency_supportlabbench.host.LogStderr.settings attribute, 52](#)
[ConcurrentException, 59](#)
[concurrently\(\)in module labbench.util, 57](#)
[ConcurrentRunnerclass in labbench.util, 59](#)
[ConfigStoreclass in labbench.util, 59](#)
[configure\(\)labbench.backends.SerialLoggingDevice method, 15](#)
[connect\(\)labbench.backends.CommandLineWrapper method, 4](#)
[connect\(\)labbench.backends.DotNetDevice method, 6](#)
[connect\(\)labbench.backends.EmulatedVISADevice method, 8](#)
[connect\(\)labbench.backends.LabviewSocketInterface method, 10](#)
[connect\(\)labbench.backends.SerialDevice method, 12](#)
[connect\(\)labbench.backends.SerialLoggingDevice method, 15](#)
[connect\(\)labbench.backends.TelnetDevice method, 18](#)
[connect\(\)labbench.backends.VISADevice method, 20](#)
[connect\(\)labbench.backends.Win32ComDevice method, 24](#)
[connect\(\)labbench.core.Device method, 31](#)
[connect\(\)labbench.host.Email method, 50](#)
[connect\(\)labbench.host.Host method, 47](#)
[connect\(\)labbench.host.LogStderr method, 52](#)
[connectedlabbench.backends.CommandLineWrapper.state attribute, 6](#)
[connectedlabbench.backends.DotNetDevice.state attribute, 7](#)
[connectedlabbench.backends.EmulatedVISADevice.state attribute, 9](#)
[connectedlabbench.backends.LabviewSocketInterface.state attribute, 11](#)
[connectedlabbench.backends.SerialDevice.state attribute, 14](#)
[connectedlabbench.backends.SerialLoggingDevice.state attribute, 18](#)
[connectedlabbench.backends.TelnetDevice.state attribute, 20](#)
[connectedlabbench.backends.VISADevice.state attribute, 23](#)
[connectedlabbench.backends.Win32ComDevice.state attribute, 25](#)
[connectedlabbench.core.Device.state attribute, 32](#)
[connectedlabbench.host.Email.state attribute, 51](#)
[connectedlabbench.host.Host.state attribute, 49](#)
[connectedlabbench.host.LogStderr.state attribute, 53](#)
[ConnectionError, 27](#)
[copy\(\)labbench.util.FilenameDict method, 59](#)
[data_formatlabbench.backends.SerialLoggingDevice.settings attribute, 16](#)
[default_valuelabbench.core.Bool attribute, 29](#)
[default_valuelabbench.core.List attribute, 29](#)
[default_valuelabbench.core.Unicode attribute, 28](#)
[default_value_repr\(\)labbench.core.Dict method, 30](#)
[default_value_repr\(\)labbench.core.List method, 29](#)
[define\(\)labbench.backends.CommandLineWrapper.settings class method, 5](#)
[define\(\)labbench.backends.DotNetDevice.settings class method, 7](#)
[define\(\)labbench.backends.EmulatedVISADevice.settings class method, 9](#)
[define\(\)labbench.backends.LabviewSocketInterface.settings class method, 11](#)
[define\(\)labbench.backends.SerialDevice.settings class method, 13](#)
[define\(\)labbench.backends.SerialLoggingDevice.settings class method, 16](#)
[define\(\)labbench.backends.TelnetDevice.settings class method, 19](#)
[define\(\)labbench.backends.VISADevice.settings class method, 22](#)
[define\(\)labbench.backends.Win32ComDevice.settings class method, 25](#)
[define\(\)labbench.core.Device.settings class method, 32](#)
[define\(\)labbench.host.Email.settings class method, 50](#)
[define\(\)labbench.host.Host.settings class method, 48](#)
[define\(\)labbench.host.LogStderr.settings class method, 52](#)
[delaylabbench.backends.LabviewSocketInterface attribute, 10](#)
[devalidate\(\)labbench.core.CaselessBytesEnum method, 29](#)

Deviceclass in labbench.core, 30
Device.settingsclass in labbench.core, 31
Device.stateclass in labbench.core, 32
DeviceConnectionLost, 27
DeviceException, 27
DeviceFatalError, 27
DeviceNotReady, 27
deviceslabbench.notebooks.panel attribute, 55
DeviceStateError, 28
Dictclass in labbench.core, 29
disconnect()labbench.backends.CommandLineWrapper method, 4
disconnect()labbench.backends.DotNetDevice method, 6
disconnect()labbench.backends.EmulatedVISADevice method, 8
disconnect()labbench.backends.LabviewSocketInterface method, 10
disconnect()labbench.backends.SerialDevice method, 12
disconnect()labbench.backends.SerialLoggingDevice method, 15
disconnect()labbench.backends.TelnetDevice method, 18
disconnect()labbench.backends.VISADevice method, 21
disconnect()labbench.backends.Win32ComDevice method, 24
disconnect()labbench.core.Device method, 31
disconnect()labbench.host.Email method, 50
disconnect()labbench.host.Host method, 47
disconnect()labbench.host.LogStderr method, 52
dll_namelabbench.backends.DotNetDevice attribute, 6
doc_attrslabbench.core.CaselessBytesEnum attribute, 29
doc_attrslabbench.core.CaselessStrEnum attribute, 30
doc_attrslabbench.core.Float attribute, 28
doc_attrslabbench.core.Int attribute, 28
DotNetDeviceclass in labbench.backends, 6
DotNetDevice.settingsclass in labbench.backends, 7
DotNetDevice.stateclass in labbench.backends, 7
dsrdtrlabbench.backends.SerialDevice.settings attribute, 14
dsrdtrlabbench.backends.SerialLoggingDevice.settings attribute, 17

element_error()labbench.core.Dict method, 30
element_error()labbench.core.List method, 29
Emailclass in labbench.host, 49
Email.settingsclass in labbench.host, 50
Email.stateclass in labbench.host, 51
EmulatedVISADeviceclass in labbench.backends, 8
EmulatedVISADevice.settingsclass in labbench.backends, 8
EmulatedVISADevice.stateclass in labbench.backends, 9
error()labbench.core.Dict method, 30
error()labbench.core.List method, 29
exception_on_stderrlabbench.backends.CommandLineWrapper attribute, 4
failure_messagelabbench.host.Email.settings attribute, 50
fetch()labbench.backends.SerialLoggingDevice method, 15
FilenameDictclass in labbench.util, 59
Floatclass in labbench.core, 28
flush()labbench.host.LogStderr method, 52
foreground()labbench.backends.CommandLineWrapper method, 4
frame()labbench.util.ConfigStore class method, 59
from_filename()labbench.util.FilenameDict class method, 60
from_hwid()labbench.backends.SerialDevice class method, 13
from_hwid()labbench.backends.SerialLoggingDevice class method, 15
from_index()labbench.util.FilenameDict class method, 60
fromkeys()labbench.util.FilenameDict class method, 60

generatorslabbench.backends.EmulatedVISADevice attribute, 8
get()labbench.data.StateAggregator method, 35
get()labbench.data.StatesToCSV method, 39
get()labbench.data.StatesToRelationalTable method, 37
get()labbench.data.StatesToSQLite method, 42
get()labbench.util.FilenameDict method, 60
get_git_browse_urllabbench.host.Host attribute, 47
get_git_commitlabbench.host.Host attribute, 47
get_git_remote_urllabbench.host.Host attribute, 47
get_hostnamelabbench.host.Host attribute, 47
get_loglabbench.host.Host attribute, 47
get_timelabbench.host.Host attribute, 47
getter()labbench.backends.CommandLineWrapper.state class method, 6
getter()labbench.backends.DotNetDevice.state class method, 8
getter()labbench.backends.EmulatedVISADevice.state class method, 9
getter()labbench.backends.LabviewSocketInterface.state class method, 11
getter()labbench.backends.SerialDevice.state class method, 14
getter()labbench.backends.SerialLoggingDevice.state class method, 18
getter()labbench.backends.TelnetDevice.state class method, 20
getter()labbench.backends.VISADevice.state class method, 23
getter()labbench.backends.Win32ComDevice.state class method, 25
getter()labbench.core.Device.state class method, 32
getter()labbench.host.Email.state class method, 51
getter()labbench.host.Host.state class method, 49
getter()labbench.host.LogStderr.state class method, 53
git_browse_urllabbench.host.Host.state attribute, 49

[git_commit_idlabbench.host.Host.state attribute, 49](#)
[git_remote_urllabbench.host.Host.state attribute, 49](#)

[hash_caller\(\)in module labbench.util, 62](#)
[Hostclass in labbench.host, 47](#)
[Host.settingsclass in labbench.host, 47](#)
[Host.stateclass in labbench.host, 48](#)
[hostnamelabbench.host.Host.state attribute, 49](#)

[identitylabbench.backends.EmulatedVISADevice.state attribute, 9](#)
[identitylabbench.backends.VISADevice.state attribute, 23](#)
[iloclabbench.util.FilenameDict attribute, 60](#)
[index_labellabbench.data.StatesToCSV attribute, 39](#)
[index_labellabbench.data.StatesToRelationalTable attribute, 37](#)
[index_labellabbench.data.StatesToSQLite attribute, 42](#)
[info\(\)labbench.core.CaselessStrEnum method, 30](#)
[info\(\)labbench.core.Dict method, 30](#)
[info\(\)labbench.core.List method, 29](#)
[info_textlabbench.core.Bool attribute, 29](#)
[info_textlabbench.core.Bytes attribute, 28](#)
[info_textlabbench.core.CaselessBytesEnum attribute, 29](#)
[info_textlabbench.core.Complex attribute, 28](#)
[info_textlabbench.core.Float attribute, 28](#)
[info_textlabbench.core.Int attribute, 28](#)
[info_textlabbench.core.TCPAddress attribute, 30](#)
[info_textlabbench.core.Unicode attribute, 28](#)
[instance_init\(\)labbench.core.Dict method, 30](#)
[instance_init\(\)labbench.core.List method, 29](#)
[Intclass in labbench.core, 28](#)
[items\(\)labbench.util.FilenameDict method, 60](#)

[keylabbench.util.FilenameDict attribute, 60](#)
[key\(\)labbench.data.StateAggregator method, 35](#)
[key\(\)labbench.data.StatesToCSV method, 39](#)
[key\(\)labbench.data.StatesToRelationalTable method, 37](#)
[key\(\)labbench.data.StatesToSQLite method, 42](#)
[keys\(\)labbench.util.FilenameDict method, 60](#)
[kill\(\)labbench.backends.CommandLineWrapper method, 4](#)
[kill_by_name\(\)in module labbench.util, 62](#)
[klasslabbench.core.Dict attribute, 30](#)
[klasslabbench.core.List attribute, 29](#)

[labbenchmodule, 69](#)
[labbench.backendsmodule, 3](#)
[labbench.coremodule, 27](#)
[labbench.datamodule, 35](#)
[labbench.hostmodule, 47](#)
[labbench.notebooksmodule, 55](#)
[labbench.utilmodule, 57](#)
[labbench.versionmodule, 67](#)
[LabviewSocketInterfaceclass in labbench.backends, 10](#)

[LabviewSocketInterface.settingsclass in labbench.backends, 10](#)
[LabviewSocketInterface.stateclass in labbench.backends, 11](#)
[length_error\(\)labbench.core.List method, 29](#)
[librarylabbench.backends.DotNetDevice attribute, 6](#)
[linspace\(\)in module labbench.notebooks, 56](#)
[Listclass in labbench.core, 29](#)
[list_devices\(\)in module labbench.core, 33](#)
[list_ports\(\)labbench.backends.SerialDevice static method, 13](#)
[list_ports\(\)labbench.backends.SerialLoggingDevice static method, 16](#)
[list_resources\(\)labbench.backends.VISADevice class method, 21](#)

[loglabbench.host.Host.state attribute, 49](#)
[loglabbench.host.LogStderr attribute, 52](#)
[log_formatlabbench.host.Host attribute, 47](#)
[log_progress\(\)in module labbench.notebooks, 55](#)
[LogStderrclass in labbench.host, 51](#)
[LogStderr.settingsclass in labbench.host, 52](#)
[LogStderr.stateclass in labbench.host, 52](#)

[make\(\)labbench.util.Testbed method, 64](#)
[make_dynamic_default\(\)labbench.core.Dict method, 30](#)
[make_dynamic_default\(\)labbench.core.List method, 29](#)
[master_filenamelabbench.data.StatesToSQLite attribute, 42](#)

[max_queue_sizelabbench.backends.SerialLoggingDevice.settings attribute, 17](#)
[metadata\(\)labbench.host.Host method, 47](#)

[ncolslabbench.notebooks.panel attribute, 55](#)
[no_state_argumentslabbench.backends.CommandLineWrapper attribute, 4](#)
[nonscalar_file_typelabbench.data.StatesToCSV attribute, 40](#)

[observe\(\)labbench.data.StateAggregator method, 35](#)
[observe\(\)labbench.data.StatesToCSV method, 40](#)
[observe\(\)labbench.data.StatesToRelationalTable method, 37](#)
[observe\(\)labbench.data.StatesToSQLite method, 42](#)
[observe_settings\(\)labbench.data.StateAggregator method, 35](#)
[observe_settings\(\)labbench.data.StatesToCSV method, 40](#)
[observe_settings\(\)labbench.data.StatesToRelationalTable method, 37](#)
[observe_settings\(\)labbench.data.StatesToSQLite method, 42](#)
[observe_states\(\)labbench.data.StateAggregator method, 35](#)
[observe_states\(\)labbench.data.StatesToCSV method, 40](#)

`observe_states()`labbench.data.StatesToRelationalTable method, 37

`observe_states()`labbench.data.StatesToSQLite method, 42

`open()`labbench.data.StatesToCSV method, 40

`open()`labbench.data.StatesToRelationalTable method, 38

`open()`labbench.data.StatesToSQLite method, 43

`optionslabbench.backends.EmulatedVISADevice.state` attribute, 10

`optionslabbench.backends.VISADevice.state` attribute, 23

`overlap_and_block()`labbench.backends.VISADevice method, 21

`panelclass` in labbench.notebooks, 55

`paritylabbench.backends.SerialDevice.settings` attribute, 14

`paritylabbench.backends.SerialLoggingDevice.settings` attribute, 17

`peekitem()`labbench.util.FilenameDict method, 60

`poll_ratelabbench.backends.SerialLoggingDevice.settings` attribute, 17

`pop()`labbench.util.FilenameDict method, 61

`popitem()`labbench.util.FilenameDict method, 61

`portlabbench.backends.TelnetDevice.settings` attribute, 19

`preset()`labbench.backends.VISADevice method, 21

`query()`labbench.backends.VISADevice method, 21

`range()`in module labbench.notebooks, 55

`read()`in module labbench.data, 44

`read()`labbench.backends.LabviewSocketInterface method, 10

`read_relational()`in module labbench.data, 44

`read_stdout()`labbench.backends.CommandLineWrapper method, 4

`read_terminationlabbench.backends.VISADevice.settings` attribute, 22

`recipientslabbench.host.Email.settings` attribute, 50

`resourcelabbench.backends.CommandLineWrapper.settings` attribute, 5

`resourcelabbench.backends.DotNetDevice.settings` attribute, 7

`resourcelabbench.backends.EmulatedVISADevice.settings` attribute, 9

`resourcelabbench.backends.LabviewSocketInterface` attribute, 10

`resourcelabbench.backends.LabviewSocketInterface.settings` attribute, 11

`resourcelabbench.backends.SerialDevice.settings` attribute, 14

`resourcelabbench.backends.SerialLoggingDevice.settings` attribute, 17

`resourcelabbench.backends.TelnetDevice.settings` attribute, 19

`resourcelabbench.backends.VISADevice.settings` attribute, 22

`resourcelabbench.backends.Win32ComDevice.settings` attribute, 25

`resourcelabbench.core.Device.settings` attribute, 32

`resourcelabbench.host.Email.settings` attribute, 51

`resourcelabbench.host.Host.settings` attribute, 48

`resourcelabbench.host.LogStderr.settings` attribute, 52

`respawnlabbench.backends.CommandLineWrapper` attribute, 4

`retry()`in module labbench.util, 63

`rtsctslabbench.backends.SerialDevice.settings` attribute, 14

`rtsctslabbench.backends.SerialLoggingDevice.settings` attribute, 17

`running()`labbench.backends.CommandLineWrapper method, 4

`running()`labbench.backends.SerialLoggingDevice method, 16

`rx_buffer_sizelabbench.backends.LabviewSocketInterface` attribute, 10

`rx_portlabbench.backends.LabviewSocketInterface` attribute, 10

`send_summary()`labbench.host.Email method, 50

`senderlabbench.host.Email.settings` attribute, 51

`sequentially()`in module labbench.util, 58

`SerialDeviceclass` in labbench.backends, 12

`SerialDevice.settingsclass` in labbench.backends, 13

`SerialDevice.stateclass` in labbench.backends, 14

`SerialLoggingDeviceclass` in labbench.backends, 15

`SerialLoggingDevice.settingsclass` in labbench.backends, 16

`SerialLoggingDevice.stateclass` in labbench.backends, 17

`set_backend()`labbench.backends.EmulatedVISADevice class method, 8

`set_backend()`labbench.backends.VISADevice class method, 21

`set_device_labels()`labbench.data.StateAggregator method, 36

`set_device_labels()`labbench.data.StatesToCSV method, 40

`set_device_labels()`labbench.data.StatesToRelationalTable method, 38

`set_device_labels()`labbench.data.StatesToSQLite method, 43

`set_path_format()`labbench.data.StatesToCSV method, 41

`set_path_format()`labbench.data.StatesToRelationalTable method, 38

`set_path_format()`labbench.data.StatesToSQLite method, 43

`set_queue()`labbench.util.Call method, 59

`set_relational_file_format()`labbench.data.StatesToCSV method, 41

[set_relational_file_format\(\)](#)labbench.data.StatesToRelationalTable method, 38
[set_relational_file_format\(\)](#)labbench.data.StatesToSQLite method, 43
[set_row_preprocessor\(\)](#)labbench.data.StatesToCSV method, 41
[set_row_preprocessor\(\)](#)labbench.data.StatesToRelationalTable method, 38
[set_row_preprocessor\(\)](#)labbench.data.StatesToSQLite method, 43
[setdefault\(\)](#)labbench.util.FilenameDict method, 62
[setter\(\)](#)labbench.backends.CommandLineWrapper.state class method, 6
[setter\(\)](#)labbench.backends.DotNetDevice.state class method, 8
[setter\(\)](#)labbench.backends.EmulatedVISADevice.state class method, 10
[setter\(\)](#)labbench.backends.LabviewSocketInterface.state class method, 12
[setter\(\)](#)labbench.backends.SerialDevice.state class method, 14
[setter\(\)](#)labbench.backends.SerialLoggingDevice.state class method, 18
[setter\(\)](#)labbench.backends.TelnetDevice.state class method, 20
[setter\(\)](#)labbench.backends.VISADevice.state class method, 23
[setter\(\)](#)labbench.backends.Win32ComDevice.state class method, 26
[setter\(\)](#)labbench.core.Device.state class method, 32
[setter\(\)](#)labbench.host.Email.state class method, 51
[setter\(\)](#)labbench.host.Host.state class method, 49
[setter\(\)](#)labbench.host.LogStderr.state class method, 53
[setup\(\)](#)labbench.data.StatesToCSV method, 41
[setup\(\)](#)labbench.data.StatesToRelationalTable method, 38
[setup\(\)](#)labbench.data.StatesToSQLite method, 44
[setup\(\)](#)labbench.util.Call static method, 59
[show_messages\(\)](#)in module labbench.util, 63
[sleep\(\)](#)in module labbench.util, 63
[start\(\)](#)labbench.backends.SerialLoggingDevice method, 17
[startup\(\)](#)labbench.util.Testbed method, 64
[StateAggregator](#)class in labbench.data, 35
[StatesToCSV](#)class in labbench.data, 39
[StatesToRelationalTable](#)class in labbench.data, 36
[StatesToSQLite](#)class in labbench.data, 41
[status_bytelabbench.backends.EmulatedVISADevice.state attribute, 10
\[status_bytelabbench.backends.VISADevice.state attribute, 23
\\[stop\\\(\\\)\\]\\(#\\)labbench.backends.SerialLoggingDevice method, 18
\\[stop_timeout\\]\\(#\\)labbench.backends.SerialLoggingDevice.settings attribute, 17
\\[stopbits\\]\\(#\\)labbench.backends.SerialLoggingDevice.settings attribute, 17
\\[stopwatch\\\(\\\)\\]\\(#\\)in module labbench.util, 63
\\[success_message\\]\\(#\\)labbench.host.Email.settings attribute, 51
\\[table_name\\]\\(#\\)labbench.data.StatesToSQLite attribute, 44
\\[TCPAddress\\]\\(#\\)class in labbench.core, 30
\\[TelnetDevice\\]\\(#\\)class in labbench.backends, 18
\\[TelnetDevice.settings\\]\\(#\\)class in labbench.backends, 19
\\[TelnetDevice.state\\]\\(#\\)class in labbench.backends, 19
\\[Testbed\\]\\(#\\)class in labbench.util, 64
\\[ThreadEndedByMaster\\]\\(#\\), 64
\\[ThreadSandbox\\]\\(#\\)class in labbench.util, 64
\\[timelabbench.host.Host.state attribute, 49
\\\[time_format\\\]\\\(#\\\)labbench.host.Host attribute, 49
\\\[timeout\\\]\\\(#\\\)labbench.backends.CommandLineWrapper.settings attribute, 5
\\\[timeout\\\]\\\(#\\\)labbench.backends.LabviewSocketInterface attribute, 12
\\\[timeout\\\]\\\(#\\\)labbench.backends.SerialDevice.settings attribute, 14
\\\[timeout\\\]\\\(#\\\)labbench.backends.SerialLoggingDevice.settings attribute, 17
\\\[timeout\\\]\\\(#\\\)labbench.backends.TelnetDevice.settings attribute, 19
\\\[to_feather\\\\(\\\\)\\\]\\\(#\\\)in module labbench.data, 45
\\\[tx_port\\\]\\\(#\\\)labbench.backends.LabviewSocketInterface attribute, 12
\\\[Unicode\\\]\\\(#\\\)class in labbench.core, 28
\\\[until_timeout\\\\(\\\\)\\\]\\\(#\\\)in module labbench.util, 65
\\\[update\\\\(\\\\)\\\]\\\(#\\\)labbench.util.FilenameDict method, 62
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Bool method, 29
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Bytes method, 28
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.CaselessBytesEnum method, 29
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.CaselessStrEnum method, 30
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Complex method, 28
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Dict method, 30
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Float method, 28
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Int method, 28
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.List method, 29
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.TCPAddress method, 30
\\\[validate\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Unicode method, 28
\\\[validate_elements\\\\(\\\\)\\\]\\\(#\\\)labbench.core.Dict method, 30
\\\[validate_elements\\\\(\\\\)\\\]\\\(#\\\)labbench.core.List method, 29
\\\[values\\\\(\\\\)\\\]\\\(#\\\)labbench.util.FilenameDict method, 62
\\\[VISADevice\\\]\\\(#\\\)class in labbench.backends, 20
\\\[VISADevice.settings\\\]\\\(#\\\)class in labbench.backends, 22
\\\[VISADevice.state\\\]\\\(#\\\)class in labbench.backends, 22
\\\[VISADevice.suppress_timeout\\\]\\\(#\\\)class in labbench.backends, 23\\]\\(#\\)\]\(#\)](#)

VISADevice.VISADeviceclass in labbench.backends, [20](#)

wait()labbench.backends.VISADevice method, [24](#)

widgetlabbench.notebooks.panel attribute, [55](#)

Win32ComDeviceclass in labbench.backends, [24](#)

Win32ComDevice.settingsclass in labbench.backends, [24](#)

Win32ComDevice.stateclass in labbench.backends, [25](#)

with_traceback()labbench.core.CommandNotImplementedError
method, [33](#)

with_traceback()labbench.core.ConnectionError method,
[27](#)

with_traceback()labbench.core.DeviceConnectionLost
method, [28](#)

with_traceback()labbench.core.DeviceException method,
[27](#)

with_traceback()labbench.core.DeviceFatalError method,
[27](#)

with_traceback()labbench.core.DeviceNotReady method,
[27](#)

with_traceback()labbench.core.DeviceStateError method,
[28](#)

with_traceback()labbench.util.ConcurrentException
method, [59](#)

with_traceback()labbench.util.ThreadEndedByMaster
method, [65](#)

write()labbench.backends.LabviewSocketInterface
method, [12](#)

write()labbench.backends.VISADevice method, [24](#)

write()labbench.data.StatesToCSV method, [41](#)

write()labbench.data.StatesToRelationalTable method, [38](#)

write()labbench.data.StatesToSQLite method, [44](#)

write()labbench.host.LogStderr method, [53](#)

write_stdin()labbench.backends.CommandLineWrapper
method, [6](#)

write_terminationlabbench.backends.SerialDevice.settings
attribute, [14](#)

write_terminationlabbench.backends.SerialLoggingDevice.settings
attribute, [17](#)

write_terminationlabbench.backends.VISADevice.settings
attribute, [22](#)

xonxofflabbench.backends.SerialDevice.settings attribute,
[14](#)

xonxofflabbench.backends.SerialLoggingDevice.settings
attribute, [17](#)